



### Product description

LoFric® Hydro-Kit™, single use, hydrophilic urinary catheter with salt solution for activation, intended for intermittent urinary catheterisation.

### Sterilisation and shelf life

Sterilised by irradiation. Use before expiry date on package.

### Package

30 catheters per customer box. 2 customer boxes per shipping box.

### Material composition

**Hydrophilic surface coating:** Polyvinyl pyrrolidone (PVP).

**Catheter:** Polyolefin-based elastomer (POBE).

**Sachet:** Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

**Item packaging:** Polyethylene butylacrylate (EBA).

### Product range

Reference	Type	Length	Size
421 08 253	LoFric® Hydro-Kit™ Nelaton (paediatric)	20 cm	CH 08
421 10 253	LoFric® Hydro-Kit™ Nelaton (paediatric)	20 cm	CH 10
423 08 253	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 08
423 10 253	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 10
423 12 253	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 12
423 14 253	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 14
423 16 253	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 16
423 18 253	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 18
420 08 253	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 08
420 10 253	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 10
420 12 253	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 12
420 14 253	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 14
420 16 253	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 16
420 18 253	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 18
425 10 253	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 10
425 12 253	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 12
425 14 253	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 14
425 16 253	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 16
425 18 253	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 18



User-friendly packaging.



**LoFric® Hydro-Kit™**

## References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimise the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcomes in a sustainable way, now and for the future.

**Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.co.uk](https://www.wellspect.co.uk).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.co.uk](https://www.wellspect.co.uk)

Wellspect Limited, Stroudwater Business Park, Brunel Way, Stonehouse, Gloucestershire, GL10 3GB, Great Britain.  
Phone: 0800 652 3350. [www.wellspect.co.uk](https://www.wellspect.co.uk)

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

CE 2797

Rx Only