



Peel-open tab.

Adhesive tab on back for reliable hanging.

Instant activation - sterile salt solution included, squeeze and it's ready to use.

Adjustable Insertion Grip allows better grip and control, promoting a non-touch technique.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{1,2}

Smooth eyelets.^{3,4}



Foldable into handy, discreet pocket size.



Nelaton Tiemann Flexible

Product description

LoFric® Origo™, single use, hydrophilic urinary catheter with salt solution for activation.

Sterilisation and shelf life

Sterilised by irradiation. Use before expiry date on package.

Package

30 catheters per customer box. 4 customer boxes per shipping box.
5 catheters per sample pack.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sachet: Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

Item packaging:

Top foil: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

Bottom foil: Polyolefin.

Product range

Reference	Type	Length	Size
442 08 25	LoFric® Origo™ Nelaton (paediatric)	30 cm	CH 08
442 10 25	LoFric® Origo™ Nelaton (paediatric)	30 cm	CH 10
442 12 25	LoFric® Origo™ Nelaton (paediatric)	30 cm	CH 12
443 10 25	LoFric® Origo™ Nelaton (male)	40 cm	CH 10
443 12 25	LoFric® Origo™ Nelaton (male)	40 cm	CH 12
443 14 25	LoFric® Origo™ Nelaton (male)	40 cm	CH 14
443 16 25	LoFric® Origo™ Nelaton (male)	40 cm	CH 16
443 18 25	LoFric® Origo™ Nelaton (male)	40 cm	CH 18
444 10 25	LoFric® Origo™ Tiemann (male)	40 cm	CH 10
444 12 25	LoFric® Origo™ Tiemann (male)	40 cm	CH 12
444 14 25	LoFric® Origo™ Tiemann (male)	40 cm	CH 14
444 16 25	LoFric® Origo™ Tiemann (male)	40 cm	CH 16
444 18 25	LoFric® Origo™ Tiemann (male)	40 cm	CH 18
445 12 25	LoFric® Origo™ Flexible (male)	40 cm	CH 12
445 14 25	LoFric® Origo™ Flexible (male)	40 cm	CH 14

All products available in sample pack, 5 catheters.
To order replace the last two digits in the Reference number with 05.



User-friendly packaging.



Sample pack, 5 catheters.

LoFric® Origo™



References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimise the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcomes in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.co.uk](https://www.wellspect.co.uk).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.co.uk](https://www.wellspect.co.uk)

Wellspect Limited, Stroudwater Business Park, Brunel Way, Stonehouse, Gloucestershire, GL10 3GB, Great Britain.
Phone: 0800 652 3350. www.wellspect.co.uk

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only