



The funnel has a different color depending on the diameter of the catheter.

Adhesive tab on back for reliable hanging.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.<sup>1,2</sup>

Smooth eyelets.<sup>3,4</sup>



Nelaton Tiemann



User-friendly packaging.

### Product description

LoFric®, single-use, hydrophilic urinary catheter, intended for intermittent urinary catheterisation.

### Sterilisation and shelf life

Sterilised by ethylene oxide gas. Use before expiry date on package.

### Package

30 catheters per customer box. 2 or 4 customer boxes per shipping box.

### Material composition

**Hydrophilic surface coating:** Polyvinyl pyrrolidone (PVP) and sodium chloride (NaCl).

**Catheter:** Polyolefin-based elastomer (POBE).

**Item packaging:** Lacquered paper, laminated film of Polypropylene (PP)/Polyamide (PA) / Polyethylene (PE).

### Product range

Reference	Type	Length	Size
400 08 25	LoFric® Nelaton (male)	40 cm	CH 08
400 10 25	LoFric® Nelaton (male)	40 cm	CH 10
400 12 25	LoFric® Nelaton (male)	40 cm	CH 12
400 14 25	LoFric® Nelaton (male)	40 cm	CH 14
400 16 25	LoFric® Nelaton (male)	40 cm	CH 16
400 18 25	LoFric® Nelaton (male)	40 cm	CH 18
400 20 25	LoFric® Nelaton (male)	40 cm	CH 20
400 22 25	LoFric® Nelaton (male)	40 cm	CH 22
401 06 25	LoFric® Nelaton (paediatric)	20 cm	CH 06
402 08 25	LoFric® Nelaton (paediatric)	30 cm	CH 08
402 10 25	LoFric® Nelaton (paediatric)	30 cm	CH 10
403 08 25	LoFric® Nelaton (female)	20 cm	CH 08
403 10 25	LoFric® Nelaton (female)	20 cm	CH 10
403 12 25	LoFric® Nelaton (female)	20 cm	CH 12
403 14 25	LoFric® Nelaton (female)	20 cm	CH 14
403 16 25	LoFric® Nelaton (female)	20 cm	CH 16
403 18 25	LoFric® Nelaton (female)	20 cm	CH 18
404 08 25	LoFric® Nelaton (female)	15 cm	CH 08
404 10 25	LoFric® Nelaton (female)	15 cm	CH 10
404 12 25	LoFric® Nelaton (female)	15 cm	CH 12
404 14 25	LoFric® Nelaton (female)	15 cm	CH 14
405 10 25	LoFric® Tiemann (male)	40 cm	CH 10
405 12 25	LoFric® Tiemann (male)	40 cm	CH 12
405 14 25	LoFric® Tiemann (male)	40 cm	CH 14
405 16 25	LoFric® Tiemann (male)	40 cm	CH 16
405 18 25	LoFric® Tiemann (male)	40 cm	CH 18

## References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimise the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcomes in a sustainable way, now and for the future.

**Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.co.uk](https://www.wellspect.co.uk).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.co.uk](https://www.wellspect.co.uk)

Wellspect Limited, Stroudwater Business Park, Brunel Way, Stonehouse, Gloucestershire, GL10 3GB, Great Britain.  
Phone: 0800 652 3350. [www.wellspect.co.uk](https://www.wellspect.co.uk)

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

CE 2797

Rx Only